

Product Name :
High Frequency Mobile Digital C-arm System**Product Code :**
NMCBS28011**Description :**

High Frequency Mobile Digital C-arm System Use : It is widely used in orthopedics, urology surgery, spinal surgery, abdominal, pain management, digestive system, gynecology and obstetrics etc. Features : High power generator with digital CCD camera to provide high resolution image Pulse fluoroscopy with controlled exposure and ultra-low radiation dose. Double high resolution LCD monitors ensure image quality . Work station has high-capacity digital storage function. Fluoroscopy and digital spot film can be stored in digitized form. It also has functions such as edge enhancement, multiple mirror image, gamma correction, cine loop, window width & window level, expert template and CD writer. It has dicom 3.0 interface to connect the internet Electric motor control for accurate positioning and all round rotation . Dual Focus : Small/Large : 0.3/0.6 Inverter Frequency : 60kHz

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,