

**Product Name :**  
Ventilator with Air Compressor**Product Code :**  
ICU0002**Description :**

Ventilator with Air Compressor Work condition : Gas source : O<sub>2</sub> (For medical use), Air Compressor (For medical use) Rated Working Pressure : 280kPa ~ 600kPa Power Source : Voltage : ~220V  $\pm$  1/2 22V, Frequency : 50Hz  $\pm$  1/2 1Hz Ventilation Mode : IPPV, A/C, SIMC, SPONT, CPAP, PSV

**Technical Specification :**

Artificial respiration is a high-performance system for providing Clean and Dry compressed air for a wide variety of medical devices & applications, such as ventilators, CPAP, BiPAP, etc.

Quiet Operation, suitable for bedside use

Compact & Modular Design

Oil-free compressor motor with thermal cutout

Automatic water removal system

Large locking castors for easy transport

Metallic cabinet with bumpers

Optional digital display with alarms

Optional Electronic Dryer

