

Product Name :
Bilevel Positive Airway Pressure with Multiple Mode**Product Code :**
ICU0005**Description :**

Bilevel Positive Airway Pressure with Multiple Mode Application : NET BPAP is mainly used for treating obstructive sleep apnea-hypopnea syndrome (OSAHS) and it is not for patient life support. Features: Work mode : CPAP, S, S/T, T Air leak compensation: Dynamic monitoring of the pressure compensation/ Dynamic monitoring of flow compensation. Large color LCD display: It is very convenient to see the setting parameters and statistical data. It can display waveform of breathing air and pressure in real time. Equipped with Multi-parameter Alarm, Mask off alarm, suffocation alarm, high respiratory rate alarm, low respiratory rate alarm

Technical Specification :

With Ti Min/ Ti Max RESmart GII BPAP Series provide non- invasive ventilation for patients with respiratory Insufficiency in hospital or at home.

Ti Min/ Ti Max ensure inspiration time and improve the treatment effect significantly. Real time pressure & flow wave display.

Oximeter Blood oxygen saturation and pulse rate display Eco Smart Technology Eco Smart heating system with innovative dual water chambers design.

Accurate water quantity control, with real time compensation.

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]