


Product Name :
Multi-Parameter Monitor**Product Code :**
ICU0006**Description :**

Multi-Parameter Monitor High resolution 12.1" color TFT display Lightweight, compact and portable ECG, SPO2, NIBP, RESP, 2-TEMP, PR Optional : 2-IBP, ETCO2, And thermal printer Built-in rechargeable lithium battery ECG waveforms of 7-leads display on the same screen 72-hours graphic and tabular trends of all parameters 72 alarm events of all parameters recall 32 seconds full-disclosure waveform review 500 NIBP measurement data storage and recall Data and waveforms color be adjustable Arrhythmia analysis and S-T segment analysis Suitable for adult paediatric and neonatal patient Standard: ECG, Respiration, Temperature, SPO2, NIBP, Pulse rate Optional: EtCO2, IBP

Technical Specification :

10" color TFT LCD screen, High brightness, Wide viewing, 8 waveforms can be displayed on the same screen. lightness and volume are adjustable independently.

Streamlined outline design with sense of the times & foldable handle.

Multi-language interface selection: Chinese, English, Spanish, Portuguese, French, Arabic, Russian, Italian, Turkish, etc.

User-friendly and man-machine dialogue interface with pop-up menus, dialogue boxes, rotary mouse and etc.

Low-power consumption design with standby mode, rechargeable high-energy built-in battery.

blood pressure measurement fast and accurately, with venipuncture function.

Advanced digital technology, accurately measure SpO2 in low perfusion situation.

Drug calculation, ST-segment detection, pace maker and arrhythmia analysis.

Display format include Oxy CRG diagram, 7-lead ECG waveform, big-font mode simultaneously ; can store 5 kinds of user-defined display format;

Sound and light integrated alarm, Alarm parameters can be adjusted to upper or lower.

Suitable for adults, pediatrics and neonates

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE,](#)