

**Product Name :**  
Electrosurgical Unit : 350W**Product Code :**  
OT Equipment0013**Description :**

Electrosurgical Unit : 350W Standard Matching Accessories: Disposable Electrosurgical Pencil : 5 pcs  
Disposable ESU Pad (split) : 5 pcs Disposable ESU Pad (solid) : 5 pcs Cable for Reusable Bipolar Forceps : 1  
pc Dual Pedal Switch : 1 pc Individual Pedal Switch : 1 pc Instructions Manual : 1 pc Features : Automatic  
Detection Function for Startup and Error Function Memory Function, 10 memory settings, previous power setting  
digits available after restart. Neutral electrode contact quality inspection system. Alarm is raised once it detect  
that the contact area is getting small to the certain dangerous situation. Large vessel sealing system, With argon  
system. With excellent compatibility, it connects with Cystoscopy, arthroscopy, laparoscopy and hysteroscopy  
and providing them the power source to cutting and coagulation. It also operated at the urethral resection of the  
prostate (TURP), the Arthroscopic surgery and the usage of gasification electrode surgery (TUVP) etc. Each  
output terminal is independent output design Transportation Data: Size : 43cm ( L ) x 38cm ( W ) x 16.5cm ( H )  
(UNIT) ; 56 ( L ) x 63 ( W ) x 28cm ( H ) (outer carton) Weight : 7kg (Net Weight) ; 10.5kg (Gross Weight)