

Product Name :
Pen-Type Digital Thermometer

Product Code :
NMDIG10022



Description :

Display range: 32C-42C(90F - 107.6F), Accuracy: 0.1 C(0.2F), Min Scale: 0.1, Memory: Last measuring reading, Measurement time: 60± seconds in oral, 100 seconds underarm, Beep function, Auto shut-off, Battery: 1.5V button battery (LR/SR-41). Pen-Type Digital Thermometer, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers