

**Product Name :**  
Ceiling Shadowless OT Light Dia 465mm x 2 domes**Product Code :**  
OT Equipment0031**Description :**

Ceiling Shadowless OT Light Dia 465mm x 2 domes Key Structure The Compact light head with shadow reduction control makes this light suitable for wide variety of procedures including the General Surgery. Every light head is equipped with two bulbs (the primary bulb and the back-up bulb). When the primary bulb fails, the back-up bulb lights up automatically to ensure continuous uninterrupted lighting throughout the case. Lights arms are made of Stainless Steel Fitting. During long time operation, there is no fatigue to the eyes, due to the range of the color temp - 4,200 K  $\pm$  300 as natural color tone. Sterilizable handle is provided in the middle of the dome in feather touch moment. Digital Control Panel to increase or decrease the intensity of the light. Technical specifications Dome head : 465 mm x 2 Halogen Bulbs : 24V x 150W x 4 nos. Light Intensity at 1 mt : 1,60,000 Lux Height Adjustment : 1065mm Dichoric Glass Reflector : 365mm Color Temperature : 4200 K  $\pm$  300 Temp rise in field : 2<sup>o</sup> - 3<sup>o</sup> C Control Unit : SMPS

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](#)