

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infrared Ear Thermometer

Product Code :
NMDIG10023



Description :

Infrared Ear Thermometer: Display range: 34C~44.0C(93.2F~111.2F), Accuracy: 0.2C(0.4F), Min Scale: 0.1, Memory: 10 memories, Measuring time: One Second, Beep function, Auto shut-off in 1 min, Battery: AAAx2, Alkaline batteries, Battery life: Approx 4000 times, Safe and hygienic Infrared Ear Thermometer, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)