

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mobile Shadowless OT Light with Single Dome Dia 405mm

Product Code :
OT Equipment0049



Description :

Mobile Shadowless OT Light with Single Dome Dia 405mm Key Structure The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Light arms are made of Stainless steel Fitting. Sterilizable handle has been provided in the middle of the dome. Inbuilt Parabolic Reflector, which increases the Lux up to 60,000 (10%). Heat absorbing & color correcting filter glass. Technical specifications
Dome head : 405mm x 1 Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 60,000 Lux Height Adjustment : 1014mm Parabolic Reflector : 355mm Color Temperature : 4200 K 300 Temp. rise in field : 2- 3 C Control Unit : Transformer

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)