

**Product Name :**  
Trial Frame

**Product Code :**  
OPPTH2010002



**Description :**

Hospital Trial Frame - 20 MM Reduced aperture. Hold 4 + 48 Lenses, I.P.D. Adjustment 50 mm to 78 mm Trial Frame, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers

**Technical Specification :**

Perfect for basic refraction needs, testing children or a practice back-up.  
PD adjustable using sliding bridge mechanism.  
The most simple and lightweight trial frame.  
Light Weight And Fully Compact.  
Holds 8 trial lenses.