

**Product Name :**  
Electronic Testing and Measuring Equipments**Product Code :**  
NMCBS102049**Description :**

Electronic Testing and Measuring Equipments We make amazing quality Electronic testing and measuring equipments which are known for their high accuracy and precision. the equipments are made using error free components. We provide them at market leading rates. We are manufacture, supplying & Export a wide range of Engineering Electronic Instruments like Electronic instruments, Characteristic of Vacuum Tubes (Valves) Characteristic of Semiconductor Devices, Lase Experiments, Nuclear Physics(G.M. Counter), Michelson Interferometers, Verification of Network Theorems, Verification of Taining Kit of Physics lab various Laws & Conversion of Signal, Characteristic of Rectifier, Filters Circuit & Regulators, Amplifiers Training Kit, Characteristic & Applications of Operational Amplifiers, Experiments of calculate e/m, Oscillators, Study of Multivibrators & applications of IC 555, Training Kit on Digital Electronic Experiements, Training Kit of Communication Lab, Function Generators, Digital/Analog testing & Measuring Instruments, Microwave Training Kit, Demonstration Training Kit, Microprocessor Trainer, Moving Coil Meters, ballistic Galvanometer, Resistance Box, Post Office Box, Condenser Boxes, UV VIS Spectrophotometer, Microprocessor UV-VIS Spectrophotometer, Digital Spectrophotometer, Digital Flame Photometer, Auto Karl Fischer Titrimeter, Digital PH, Conductivity, & Temperature Meters, Digital Colony Counter, Filter Photo Calorimeter, Digital Ph Meters, Poratable Water & Soil Anaalysis kits, laboratory Conductivity Meters, Total Dissolved Solids Meters, Digital Dissolved Oxygen Meters, Digital Potentiometers, Digital Turbidity Meters, Digital Microprocessor Haemoglobin Meters etc

---

Website: [SITE](#)

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**