

Product Name :
Indirect Ophthalmoscope Halogen**Product Code :**
OPPTH2010009**Description :**

Halogen Binocular Indirect Ophthalmoscope is an excellent instrument for all standard applications, as well as with uncooperative pediatric patients & in cases where the papillary diameter is less than 4mm. The need for pupillary dilatation is virtually obviated when using the SPF-capability of this dual purpose clinical instrument. Indirect Ophthalmoscope Halogen, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)