

Product Name :
Laboratory Plasticwares**Product Code :**
NMCBS102085**Description :**

Laboratory Plasticwares Ensuring th quality for all our products, we offer Lab Plasticware which are constructed using best of the materials. these plasticware have excellent creationand are highly durable and reliable. This category included manufacturers, exporters, traders, distributors of Laboratory Plasticwares e.g Animal Cage, Bottles, Trays, Flask, Pipettes, Racks, Vials, Containers, Beakers, Burette, Centrifuge Tubes, Conical Measures, Connectors, Coplin Jar, Cryo Vial, Cryo Box, Cuvettes, Desiccator, Fisher Clamp, Flask Stand, Forceps, Funnel Holder, Funnels, Gas Jar, Ice Bucket, Immunology Plates, Kipp's Apparatus, Leclanche Cell Pot, Magenta Box, Magnetic Needle Stand, MCT Box, Measuring Jugs, Measuring Cylinders, Menstrual Regulation Syringe Kit, Micro Slide Cabinet, Micro Tips & Micro Tip Box, Miscellaneous Products, Model, Model Set, Petri Dish, Pipette Box, Pipette Pump, Pipette Stand, Plantation Pots, Pneumatic Trough, P.P. Caps For Jam Bottle, Rain Gauge Cylinder, Regent Reservoir, Retort Stand, Simple Cell Pot, Slide Box, Slide Mailer, Spatula, Staining Box, Stirrer, Stop Cock, Storage Box, Test Tube Cap, Test Tube Stand & Basket, Test Tubes Transparent, Vernier Calipers, Vertex Mixer etc..