

Product Name :
Direct Halogen Spot Light**Product Code :**
OPPTH2010018**Description :**

State of the art Direct Halogen Spot Light with optimal halogen illumination 12v 50w suitable for operations in offices of Gynecologists and dentists and for diagnosis in other medical fields. The focusing lens at the instruments head allows adjusting the size of the illuminated examination field. The light intensity can be adjusted by means of the rheostat. The highly flexible special arm can be moved in all directions and remains in the requested position. The examination light is available as desk model, wall model or floor model with mobile steel stand (Iron cast base), black powder coated. The practical universal clamp Allows attaching the desk model also to beds or rails in hospitals. Direct Halogen Spot Light, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers