

Product Name :
Lab Glassware**Product Code :**
NMCBS102169**Description :**

Lab Glassware We are manufacturer & Exporter of Complete laboratory Glassware which are ideal for laboratory purposes. These Borosilicate Glass and glassware are designed with great precision and offer high accuracy rate. Advanced Technocracy Inc. are manufacturer and exporter of Laboratory glassware since 1954. Such as Beaker, Measuring Cylinder, Volumetric flask, Volumetric flask amber colour, Conical Flask, Burette with stop cork, Flat Bottom Flask. Round Bottom flask, Erlenmeyer flask, Boiling flask, Three neck flask, Four neck flask, two neck flask, Amber colour flask, Flask with stopper, Iodine flask, Kjeldhal Flask, Graduated pipette, Volumetric pipette, Serological Pipette, Pasteur Pipette, Mohr type pipette, Milk Pipette, Funnel, Separating funnel, Centrifuge tube, Oil centrifuge tube, Test tube, Test tube with stopper, Culture media tube, Liebig condenser, Coil condenser, bulb Condenser, Culture tubes, Distillation Apparatus, Soxhlet extraction apparatus, Dropping funnel, Distillation flask, filter flask, Flask iodine, Flask pear shape, Flask kjeldhal, Buchner funnel with sintered disc, Dean & stark, Dessicator plain, Dessicator Vacuum, Gooch Crucible, Ganong potometer, Hoffman Voltmeter, Nessler cylinder, Reagent Bottle narrow mouth, Reagent bottle wide mouth, Reagent bottle amber colour, Weighing Bottle, Specific gravity bottle, Gas wash bottle, Automatic burette, A Class Glassware with works certificate.

