

**Product Name :**  
Foldable Wheelchair with Detachable Armrest**Product Code :**  
NMWC0003**Description :**

Wheel chair folding detachable armrest (also available in fixed armrest):- 15cm castor in front with ball bearing movement & 31cm wheel at rear with brakes. This standard wheel chair specially designed for general purpose use at Hospitals/ Nursing homes / Railways / Airports / Shopping malls and other public places. Foldable Wheelchair with Detachable Armrest, Wheel Chairs Exporters, Wheel Chairs Manufacturers