

Product Name :
Wheelchair Folding (Travel) Superior**Product Code :**
NMWC0007**Description :**

Wheel chair folding detachable armrest. 20cm castors in front with ball bearing movement & 41cm wheels with grey tyre at rear with brakes, swing away and adjustable foot rest for close approach can easily be stored and carried in small car. This superior compact & lightweight chair is suitable for traveling purpose & Rough surface. Wheelchair Folding (Travel) Superior, Wheel Chairs Exporters, Wheel Chairs Manufacturers