

Product Name :
OVEN VACCUM**Product Code :**
NMCBS11004**Description :**

OVEN VACCUM Specification :- OVEN VACCUM Chamber Size : a) 250x300 12 ltrs. b) 300x300 17 ltrs. c) 300x500 40 ltrs. Ltrs Size : a) 12 ltrs. b) 17 ltrs. c) 40 ltrs. Optional Size : i) Timer 0-24 hours. ii) Imported Thermostat German, Accuracy of 0.1 \bar{A} , $^{\circ}$ C. iii) digital temp. Indicator-cum-Controller. Iv) Microprocessor PID digital temp. Indicator-cum-controller. Double walled construction. Inner round chamber made out of stainless steel, exterior made out of mild steel duly finished in white stoving enamel/powder coated paint with mat finished colour combinations. Temperature range 50 $^{\circ}$ C to 130 $^{\circ}$ C, 3 $^{\circ}$ C. The body capable of with standing vacuum of 760 mm hg. Workable on 220V AC 50 Hz single phase. Supplied without Vacuum pump OVEN VACCUM, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers