

Product Name :
WALK-IN-INCUBATOR**Product Code :**
NMCBS11008**Description :**

WALK-IN-INCUBATOR Designed for vast number of incubation applications. Can be installed even in hallway room or wherever convenient. Design feature includes enameled mild steel sheet exterior, aluminium / stainless steel, interior ventilators, internal lighting, exhaust fan etc. Suitable insulated 300 mm sq. Inspection window at the door permits inside view with -out opening the door. Complete with panel. Pre fabricated at factory assembled at site. Complete with cord and plug to work on 220 V single phase 50 Hz power supply. Temperature range from 5[°] C above ambient to 60[°] C \pm 1[°] C. Controlled by digital indicator -cum -controller. Available sizes W x H x D a) 1.2x2.3x1.2 meter b) 2.3x2.3x2.3 mete WALK-IN-INCUBATOR, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,