

**Product Name :**  
SEROLOGICAL WATER BATH (Thermostatic Water Bath)**Product Code :**  
NMCBS11009**Description :**

SEROLOGICAL WATER BATH (Thermostatic Water Bath) Temperature range from ambient temperature to 80° C. Thermostatic control with accuracy of  $\pm 0.5^\circ$  C. Double walled, inside stainless steel and outside mild steel painted in epoxy powder coating. Bath consists of two pilot lamps. Temperature control knob and on/off switch. To work on 220/230 volts A.C. Supplied with or without stirring arrangement, without racks and thermometer. Lid of water bath is made of stainless steel sheet. Size inside chamber (Stainless Steel) Capacity (Approx.) 300x200x175 mm (suitable for 2 racks) 14 ltrs. 375x300x175 mm (suitable for 4 racks) 18 ltrs. 455x300x175 mm (suitable for 6 racks) 25 ltrs. 605x300x175 mm (suitable for 8 racks) 32 ltrs. Optionals: L-shaped thermometer (fitted with water bath). Imported Thermostat "EGO" / "JUMO" German with an accuracy of  $\pm 0.5^\circ$  C Digital temp. indicator-cum-Controller. Microprocessor PID digital temp. Indicator-cum-controller. Stirrer with 1/20 H.P. High speed motor with stainless steel stirring rod and blade SEROLOGICAL WATER BATH (Thermostatic Water Bath), Heat And Refrigeration System Exporters, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers