

**Product Name :**  
OIL BATH (High temperature)**Product Code :**  
NMCBS11011**Description :**

OIL BATH (High temperature) Double wall construction. Inner chamber made of highly polished stainless steel sheet and the outer is made of mild sl duly finished in white stoving enamel/powder coated paint with mat finished colour combinations.

**Technical Specification :**

Temperature range 50½ C to 250½ C, 2½ C is controlled by a thermostat. Supplied complete with stirrer, diffuser shelves, cover etc. Specification OIL BATH Chamber Size a) 250x250x250 mm (10x10x10) b) 300x300x250 mm (12x12x10) c) 350x350x250 mm (14x14x10) Optional (i) Mercury contact thermometer and relay with accuracy better than 1½ C. (ii) Imported Thermostat with an accuracy of 1½ C (iii) Digital temp. indicator-cum-Controller. (iv) Microprocessor PID digital temp. indicator-cum-controller.