

**Product Name :**  
HEATING MANTLE**Product Code :**  
NMCBS11014**Description :**

HEATING MANTLE It with stands high temperature up to 350oC. The body of the mantle is spun in one piece from non-rusing Aluminium (up to 5.0 liters capacity) duly finished in white stoving enamel/powder coated paint with mat finished colour combinations. All heating mantels are lagged with special grade mineral wool. Works on 220/230 Volts A.C. Capacity : Wattage : 50ml 60 100ml 60 250ml 150 500ml 200 1ltr 300 2ltr 450 3ltr 450 5 ltr 600 HEATING MANTLE, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#),