

Product Name :
Melting Point Apparatus**Product Code :**
NMCBS11017**Description :**

Melting Point Apparatus For determination of melting points up to 450°C. **Manufacturers India China**

Technical Specification :

For determination of melting points up to 400°C. The mounted block is illuminated from below by a lamp fixed in its housing. The aluminium block accepts three capillary tubes and sensor. The temperature is regulated by energy regulator fitted to the unit. The lamp provides uniform and shadow less regulator fitted to the unit. No Silicon oil or other bath is required as heating media.

- Digital Temp. Indicator
- Dry Digital Melting Point
- Cooling Fan
- Temperature Hold Switch
- 3 digit LED Display
- Heat Controller: Thermostatic
- Temperature Range : 30°C - 400 o c
- Heater : 250 watts
- Power consumption: 220/230V AC Supply Single phase 50 HZ

Melting Point Apparatus Manufacturers India China, Melting Point Apparatus Lab Instruments, Medical Lab Equipment Manufacturers India China

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,