

**Product Name :**  
HORIZONTAL AUTOCLAVE**Product Code :**  
NMCBS11020**Description :**

HORIZONTAL AUTOCLAVE Most ideal for any institute or laboratory doing bulk sterilizing and research work. This autoclave has double walled chamber and a boiler and these are mounted on a heavy duty mild steel stand. The inner made of thick stainless steel sheet and Outer chamber is made of stainless steel/mild steel sheet. Automatic self closing device to prevent opening of the door while the sterilize is in operation . Stainless steel steam generation unit fitted with easily replaceable heating element. Fitted with Pressure control switch.

SPECIFICATION:- HORIZONTAL AUTOCLAVE SIZE KILLO-WATT CAPACITY 300 x 500 mm 3 KW 40 ltrs 350 x 550 mm 4 KW 50 ltrs. 400 x 600 mm 5 KW 78 ltrs. 450 x 600 mm 5 KW 98 ltrs. 500 x 750 mm 6 KW 152 ltrs. 400 x 1000 mm 9 KW 162 ltrs. PRESSURE RANGE 10 to 20 psi HORIZONTAL AUTOCLAVE, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,