

**Product Name :**  
VACCINE STORAGE CABINET**Product Code :**  
NMCBS11024**Description :**

VACCINE STORAGE CABINET Vaccine storage cabinet, the most versatile unit for vaccine for storage. Inner chamber made of stainless steel external mild steel sheet body of rigid steel frames. Insulated thoroughly provided with face fitting doors, complete with self-closing hinges lockable fasteners. Cabinet is fully sealed to avoid spilling and contamination from external environment. Temp. ranges from 4° C to 6° C. Available in various sizes for keeping vaccines bottles. Drawers with clear Perspex front and index card holders. Temperature controlled and displayed by Microprocessor PID digital temperature controlled-cum-indicator. Audio visual alarms are triggered, in the event of fluctuations exceeding permitted tolerance level. Suitable to work on 220 V, single phase, 50hz, AC supply. Available in following size:- a) 450x400x850 mm. b) 600x675x850 mm. OPTIONAL: i) Chart type Temperature recorder. ii) Voltage Stabilizer. (Recommended) VACCINE STORAGE CABINET, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,