

**Product Name :**  
DEEP FREEZER**Product Code :**  
NMCBS11025**Description :**

DEEP FREEZER Double walled with inner chamber made of S.S and outer of thick PCRC sheet duly enamel painted. Temperature range from ambient to -20° C is achieved by hermetically sealed compressor. Body is mounted on a sturdy angle iron frame and is mounted on caster wheels. Unit is fitted with solid State Digital Temperature Indicator-Cum-Controller. Horizontal Chest type models are supplied in a single compartment where as vertical models are supplied with three stainless still storing shelves. Specification DEEP FREEZER Capacity 4 cuft. 6 cuft. 10 cuft. 12 cuft. Optional Supply with automatic voltage stabilizer. Make Ajanta DEEP FREEZER, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers