

**Product Name :**  
DE-HUMIDIFIER**Product Code :**  
NMCBS11026**Description :**

DE-HUMIDIFIER, Medical Electric Portable Dehumidifier Exporters, Hospital Lab Electric Portable Dehumidifier Manufacturers

**Technical Specification :**

DE-HUMIDIFIER (Fitted with time controlled safety device in, case humidity controller fails, for safety of the compressor) These are water condensing type, based on "Refrigeration Cycle" to avoid use of expensive chemicals required for desiccant type models. The unit is self standing and mounted on a study angle iron frame. Fitted with caster wheels.

It removes excess moisture from the atmosphere and os fitted with hermetically sealed compressor. Heat exchanger and water condensing coils are made of copper. Water collection tank provided below the unit to collect condensate water.

Humidity is controlled by HUMIDITY CONTROLLER incorporated in the unit, up to 30% - 3% RH, depending upon the ambient conditions.

To achieve very low humidity, the room should have a independent window type Air conditioner to much with the dehumidifies heat load in the room. Suitable to work on 220/230 volts A.C. Supply.

**Medical Electric Portable Dehumidifier Exporters, Hospital Lab Electric Portable Dehumidifier Manufacturers**

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,