

Product Name :
Laboratory Oven**Product Code :**
NMCBS11029**Description :**

Laboratory Oven Our company is engaged in offering a wide range of laboratory ovens (Hot Air Oven) that is widely used in variety of thermal processing. These also find application in general lab work component and stability testing, core hardening, drying glassware and sterilizing. Some of the features are: • Double walled – outer mild steel & inner mild steel or stainless steel. • 80/20 Ni heating elements • Temperature range : 5°C above ambient to 300°C 2°C • 230v ac 50 Hz mains supply • Digital temperature controller Our range of laboratory oven is available in following given models: Models Chamber H x W x D (cms) Watts Shelves AEI-20a 30 x 30 x 30 750 1 AEI-20b 35 x 35 x 35 1000 2 AEI-20c 45 x 45 x 45 1500 2 AEI-20d 45 x 45 x 60 1750 2 AEI-20e 60 x 60 x 60 2000 3 AEI-20f 45 x 60 x 90 2200 3 AEI-20g 60 x 60 x 90 2500 3 Laboratory Oven, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)