

**Product Name :**  
Shaking Machines**Product Code :**  
NMCBS11030**Description :**

Shaking Machines We offer a comprehensive range of shaking machines that is commonly used in biology labs of leading educational institutions as well as research organizations. These shaking machines are platform type mounted on heavy-duty mild steel angle frame fitted with heavy cast iron pulleys. Our range of shaking machines is electronically driven by PMDC motor, which offers rotary shaking action. These have an operating speed of 200 rpm and are available with variable speed varying from 50 to 200 rpm. Some of the specifications of these shaking machines are: • Flexible spring action lotus shaped clamps • One tray is provided for holding flasks of 500 ml • Capacity: 230v ac 50hz These shaking machines are available in following models: Model Capacity To hold Flasks Platform Size L x B (cams) AEI-25a 9 x 500ml. 50 x 50 AEI-25b 12 x 500ml. 50 x 50 AEI-25c 16 x 500ml 60 x 60. AEI-25d 25 x 500ml 75 x 75 AEI-25e 36 x 500ml 80 x 80 AEI-25f 49 x 500ml 90 x 90 AEI-25g 100 x 500ml 120 x 120 Shaking Machines, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers