

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Orbital Shaking Incubators

Product Code :
NMCBS11032



Description :

Orbital Shaking Incubators Double Walled Inner SS / Outer MS; Temperature Range: 10°C to 60°C 0.5°C Preset dual display electronic digital temperature controller cum indicator; RPM Range: 50 to 250 RPM Orbital (25 cms) Variable Speed Controlled by Digital Display RPM Controller cum Indicator (Permanent Management DC Motor) Complete with St. Steel Tray with lotus shaped spring action clamps Orbital Shaking Incubators, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: SITE