

**Product Name :**  
Meat Cutting Machine**Product Code :**  
MEDICAL-ANTM10010003**Description :**

Meat Cutting Machine

**Technical Specification :**

Meat cutting machine is design for mid sized animal butcher shop and restaurants. Constructed of heavy steel angle iron clad with ms sheet duly power coated table stainless steel and used 18" blades. Electrically operated. This is ideal for those with limited space but requiring a saw capable of cutting all types of products. It has an easily removable door with smooth surface for easy cleaning and feature adjust able guide that slide down to expose only blade that is being used. Meat & Bone saw tabletop 18" s/s blade medium duty 24" cutting clearance cut. The minimum specimen will be 3mm in width. The travel of the ss table 18" in travel. Electrically operated. Stainless steel Table top. Adjustment of the cutter for upward & downward position. Total unit on casters for easy mobility.