

**Product Name :**  
Ultra Low Temperature Deep Freezer-80C**Product Code :**  
NMCBS11034**Description :**

Ultra Low Temperature Deep Freezer-80C We offer a wide range of ultra low temperature deep freezer that provides an ideal freezing environment for advanced medical and industrial applications. These are basically used in prevention of plasma, related blood components, vaccines, micro-organism, and testing of materials and electronic components at -80 Deg C. Some of the features are:

- Construction - Double Walled Construction, Exterior made of mild steel powder coated & interior made of stainless steel 304 (STD Model) or Exterior & interior made of stainless steel 304 (GMP Model)
- High density CFC Free Poly Urethane Foam insulated cabinet
- Keyed door lock for improved security with magnetic door gasket.
- Full Length inner glass door for clear inner view.
- Environment friendly Ultra Bio Freezers with micro control and CFC free cascade refrigeration system.
- Air cooled, hermetically sealed compressor unit.
- Temperature Range: -80°C controlled by microprocessor based temperature controllers.
- Temperature alarm system with adjustable high - low alarm limits includes audio visual alarm on temperature deviation.
- 230V Ac 50 HZ Mains.

Further, the technical specifications are as follows:

Order Code	Capacity Cu.ft	Ltrs	Inner Dimensions H x D x W (cms)	No. of Shelves
AEI - 5a	3 cft	100L	50 x 45 x 45	2
AEI - 5b	6 cft	180L	70 x 51 x 51	2
AEI - 5c	10 cft	300L	100 x 55 x 55	3

Ultra Low Temperature Deep Freezer-80C, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

