

**Product Name :**  
Computed Tomography CT Scan Machine**Product Code :**  
CT-SCN0002**Description :**

Computed Tomography CT Scan Machine

**Technical Specification :**

Gantry rotating faster than ever before, spectral imaging is offered relying on leading-edge hardware and software reformations. Computed Tomography CT Scan Machine is a revolutionary product.

**Features:**

Ultra-HD imaging: Micro focus spot, Quad-Sampling technology work with iHD (isotropic High Definition) and 1024 matrix, the spatial resolution can reach 30 lp/cm 0%MTF.

Uncompromised with heart rate: With 0.259s ultrafast rotation speed, Medical Computed Tomography CT Scan Machine never compromises with the high heart rate. Arrhythmia handling and CMC (Coronary Motion Artifact Clear) will also help to ensure successful coronary artery exam.

Low dose design: Unique 60kV scanning provides clinical breakthroughs on low dose scanning. It's significantly beneficial to pediatrics. Clear View iterative reconstruction algorithm can reduce dose for patients of all ages.

Dual energy innovation: Spectral imaging adds tissue characterization to morphology, improves the value of radiological department.

