

Product Name :
Medical Refrigerators**Product Code :**
MED-REF0001**Description :**

Medical Refrigerators Tender Suppliers Africa and Medical Refrigerators Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Medical Refrigerators Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Medical Refrigerators China, Medical Refrigerators Kenya.

Medical Refrigerators Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Medical Refrigerators Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),

Suriname (Paramaribo)

Technical Specification :

Temperature range/Accuracy : 2°C - 8°C , 0.1°C

Control System/ Display : Microprocessor Control, Large Screen Display

Alarm : Audio and Visual alarm for high and low temprature, door ajar, Syatem and Power failure(Except RF1050)

Refrigerant Type/Defrost : Forced Air Refrigeration System/Autodefrost

Condenser and evaporator : High Efficiency air-cooled condenser, made of Copper

Structure : Unibody Design and Mono -assemble foaming(Rigidpolyurethane insulating material)

Power Supply : 220V± 10%, 50/60 Hz

Compressor : International Famous Compressor

Size : 430x450x575

Capacity : 50L

Shelves : 1

Consumption : 50W

Refrigerant : R134A, CFC-FREE

Castors : No

Material Internal : ABS Plate

Material External : Cold-Rolled steel with anti bacteria powder coating

Door : Foam

Standard accessory : Shelves

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,