

Product Name :
Full-Functional Neonatal Nursing CPR Auscultation
Defibrillation and Pacing ECG Manikin

Product Code :
MED-CRE0002



Description :

Full-Functional Neonatal Nursing CPR Auscultation Defibrillation and Pacing ECG Manikin

Technical Specification :

Pupil observation: pupil LCD is CSTN pseudo-color, 65k color, RGB. Diameter range is 1-9mm. Pupil states: normal, mydriasis and myosis, etc;
Hair combing and face washing
Otic drops and irrigation
Oral cavity care
Airway management: realistic mouth, nose, tongue, gum, pharynx, larynx, esophagus, epiglottis, trachea and tracheal ring, can practice tracheal intubation, sputum suction, oxygen inhaling.
Oral and nasal tracheal intubation
Incision of tracheal nursing
Venipuncture, injection and transfusion in arm
Intramuscular injection in bilateral deltoid, bilateral vastus lateralis and bilateral buttock
Enema
Urethral catheterization
Holistic care: sponge washing, replacing clothes, cold and hot therapy
Limbs joints: bend, rotation, upper and lower movement
CPR operation training: support many inhaling ways, e.g.: mouth to mouth, mouth to nose, simple respirator to mouth; electronic monitoring of airway open, inhaling times, frequency and volume, and compression times, frequency, site and depth; automatically judge the ratio of artificial respiration and external chest compression; show data in real time with English voice prompts during the entire process;

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,