

**Product Name :**  
CPR Emergency Skills Training Manikin

**Product Code :**  
MED-EMG0003



**Description :**

CPR Emergency Skills Training Manikin

**Technical Specification :**

Vital signs simulation: pupil change and carotid arterial pulse

Airway open

CPR training : CPR training is designed according to 2010 international CPR guideline, can perform artificial respiration and external cardiac compression, English voice prompts during the entire process; standard airway open, real time operation waveform display, voice prompts of correct and wrong operation, data statistic ,result printing, available of training and exam mode.

Multi-media courseware teaching is available, can be connected with projector.

Available standalone version and online interactive version.

Standalone version can be upgraded to online interactive version

Standard components:

Full body manikin (1)

CPR application software (1)

CPR converter (1)

Aluminum & plastic box (1)

CPR operation pad (1)

Exchangeable lung bag (5)

Exchangeable face skin (1)

CPR face shield sheet (50pcs/box)

Temperature sensing print paper (2 roll)

Guide manual (1)



## Medical Hospital Equipments Manufacturers

**Website:**Medical Hospital Equipments Manufacturers,