

Product Name :
Half Body CPR Emergency Skills Training Manikin**Product Code :**
MED-EMG0004**Description :**

Half Body CPR Emergency Skills Training Manikin

Technical Specification :

Simulating standard open airway

External breast compression: display device and alarm device

indicator light display of correct and wrong compression ; alarm of wrong compression;

Intensity display of correct (at least 5cm) and wrong (less than 5cm) compression; alarm of wrong compression.

Artificial respiration (inhalation): display device and alarm device

Inhalation 600ml, wrong indicator light display and alarm prompting; inhalation between 500-600ml right indicator light display;

Indicator light display open airway;

Inhalation too quickly or too much result in air entering into stomach; indicator light display and alarm prompting.

Ratio of compression and artificial respiration: 30:2 (one or two person).

Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.

Operation frequency: at least 100 times per minute

Operation methods: exercise operation

Examination of pupil response: mydriasis and myosis

Examination of carotid response: pinch the pressure ball by hand and simulate carotid pulse

Working conditions: Input power is 110-240V

Medical Hospital Equipments Manufacturers

Website: [SITE](#).

NAUGRA[®]