

**Product Name :**  
Half Body CPR Emergency Skills Training Manikin

**Product Code :**  
MED-EMG0004



**Description :**

Half Body CPR Emergency Skills Training Manikin

**Technical Specification :**

Simulating standard open airway  
External breast compression: display device and alarm device  
indicator light display of correct and wrong compression ; alarm of wrong compression;  
Intensity display of correct (at least 5cm) and wrong (less than 5cm) compression; alarm of wrong compression.  
Artificial respiration (inhalation): display device and alarm device  
Inhalation 600ml, wrong indicator light display and alarm prompting; inhalation between 500-600ml right indicator light display;  
Indicator light display open airway;  
Inhalation too quickly or too much result in air entering into stomach; indicator light display and alarm prompting.  
Ratio of compression and artificial respiration: 30:2 (one or two person).  
Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.  
Operation frequency: at least 100 times per minute  
Operation methods: exercise operation  
Examination of pupil response: mydriasis and myosis  
Examination of carotid response: pinch the pressure ball by hand and simulate carotid pulse  
Working conditions: Input power is 110-240V

---

## Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**