

Product Name :
ICU Bed Electric 5 Function Column Model X-Ray
Permeable Back Rest

Product Code :
MED-ICU0002



Description :

ICU Bed Electric 5 Function Column Model X-Ray Permeable Back Rest

Technical Specification :

Electrically operated with Remote control , Back rest (0 - 77°), Foot rest (0-35°), Height adjustment, Trendelenburg Tilting (0-18°) and Reverse Trendelenburg (0-18°)
Fowler and vascular positions, Manual Knee Mechanism with Three positions
Size : L 2160mmXW955mmX H (460mm- 860 mm). Bed Extension L-180 mm.
Epoxy powder coated Framework with X-ray permeable Backrest, X-ray cassette holder .
All function controlled with Linak (Denmark) power device.
Control Panel - Remote control unit / Nurse control panel which can lock Patient remote control usage
Electronic CPR-with Nurse control panel, Manual CPR -with handles on both side of back rest (for flat position of bed backrest for emergency)
Removable head/foot boards with plastic protection bumpers (dia 90 mm) made from Polypropylene material
Foldable Side rails made from made from ABS / Polypropylene material
360° movable castors with central locking (150 mm), 5th CASTOR for easy movement of bed (125 mm).
Emergency battery backup.
Belt hook, accessories slots, urine bag holder.
Height adjustable IV pole- 2 hooks
Angle gauge to show Trendelenburg and reverse Trendelenburg angles

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,