

Product Name :
ICU Bed Super Deluxe Mechanical 5 Function**Product Code :**
MED-MEC0002**Description :**

ICU Bed Super Deluxe Mechanical 5 Function

Technical Specification :

ICU Bed (5 function) with back rest 0-75°, knee rest 0-35°, High - Low 600-800mm, Trendlenburg : 0-8° , Reverse Trendlenburg 0-8° maneuverable by 4 separate crank mechanisms from foot end.
Size: L 2030mm x W 900mm x H 600-800mm (Adjustable Height).
Epoxy coated mild steel structure with 4 section mattress platform made up of CRCA M.S. Sheet.
Detachable ABS Head and foot bows and roller bumpers.
Four easy lifting guardrails (2 on each side), which can be fixed upward and downward.
Provision for I.V. Rod on four locations.
125mm Dia. castors with diagonal locking