

**Product Name :**  
Finger Splint**Product Code :**  
MED-AMB0002**Description :**

Finger Splint

**Technical Specification :**

The product is mainly used for fastening of the 1st and 2nd joint of fingers. It is convenient to use and simple for operation. It is made of aluminium bar, sponge pad and adhesive strip. End part of the aluminium bar is bending upward, which makes the injury finger stretch upward straightly and keep fastened.

Specifications:

Color: Blue

Material: Aluminum &amp; closed cell foam

Size: S, M, L, XL