

**Product Name :**  
Urine Strip Analyzer

**Product Code :**  
MED-CLE0002



**Description :**

Urine Strip Analyzer

**Technical Specification :**

Suitable for laboratories, hospitals and clinics  
Compact, light weight and easy to use  
Audible, visual prompts  
2.8 inch color TFT display  
Three working modes for different user groups  
Compatible with 8, 10 and 11 parameter test strip  
Built-in thermal printer for realtime test results  
Dual wavelength reflectance test method  
PC communication using serial and USB port