

**Product Name :**  
Tissue Floating Bath

**Product Code :**  
MED-PTH0002



**Description :**

Tissue Floating Bath

**Technical Specification :**

Wetting temperature range: room temperature ~ 70 ° C  
Baking sheet temperature control range: room temperature ~ 100 ° C  
Control accuracy: ±2°C  
Spreading pot size: 210 x 1170 x 655 (mm)  
Baking table size: 250 x 108 (mm) (20 slides)  
Machine size: 460 x 310 x 135 (mm)