

Product Name :
Tablet Friability Tester**Product Code :**
MED-PTH0003**Description :**

Tablet Friability Tester

Technical Specification :

Two drums run synchronously. Auto stop at timing point.

The drum is made of high quality hyaline acryl.

The speed and revolution of the drum can be preset at any time, either preset or real data is displayed on the LED.

A MPU is used to control all operations, high accuracy is achieved in testing rotary speed and circle number.

Automation: auto-test, auto-diagnose, auto-alarm.

Cylinder number: 2

Cylinder radius: 286mm

Cylinder depth: 39mm

Tablet fall height: 156mm

Rotation speed: (20-90)rpm

Rotation Speed accuracy: ± 1 rpm

Rotation count: (10-900)r

Rotation count accuracy: ± 1 r

Power: 220V/ 50Hz/ 60W or 110V/ 60Hz/ 60W

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]