

**Product Name :**  
Soxhlet Extraction Fat Analyzer**Product Code :**  
MED-PTH0008**Description :**

Soxhlet Extraction Fat Analyzer

**Technical Specification :**

Small volume,fast heating speed,even heating.

Extraction bottle and apparatus are made of glass to avoid ether spills.

Accurate digital temperature controlling,easy to operate

Measuring Range: Oil content within the scope of 0.5%-60% of good,feed,oil and various kinds of fat product

Measuring Sample Quantity: 6 simultaneously

Heating time: Within 10 min.

Recovery System: Automatic

Recovery Rate: ?80%

Temperature Controlling Range: Room temperature 0-100?

Power: AC 220V/50Hz,1000W