

Product Name :
Vertical Deep Freezer**Product Code :**
NMCBS11039**Description :**

Vertical Deep Freezer We manufacture Deep Freezer for various industries Pharmaceuticals, Chemical Industries. High Quality stainless steel is used to manufacture the chambers of thick gauge and is duly polished and welded to make them completely leak proof.

Technical Specification :

These are available in two models i.e. horizontal & vertical with various capacities and temperature ranges. Some salient features of Deep Freezers manufactured by us are as follows: Cabinet in all sides 2" - 6" put insulation. Main door (Insulated with Puff) provided with silicon gasket to seal chamber from atmospheric pressure. Mild steel outer body is powder coated for lasting finish and corrosion resistance. Air cooled single stage (for -20°C and -40°C) Used with hermetically sealed compressors Controlled with digital indicator cum controller Moreover, we offer custom based freezers with Power failure alarms / temperature raise alarms/voltage stabilizer with high voltage low voltage cut-offs and other advance options. , Hospital Heat And Refrigeration Lab System Exporters, Medical Lab Heat And Refrigeration Instruments Manufacturers