

Product Name :
One Side Vertical Air Flow Clean Bench**Product Code :**
MED-PTH0011**Description :**

One Side Vertical Air Flow Clean Bench

Technical Specification :

The quasi-closed glass damper can effectively prevent the external airflow from penetrating and the odor of the human body.

The adjustable air volume fan system ensures that the working wind speed is always in an ideal state.

Economic vertical single purification workbench.

The liquid crystal panel control with LED.

The test area is made of stainless steel for ease of operation.

Overall dimensions: 850*570*1550 (W*D*H) mm

Purification area size: 700*500*500 (W*D*H) mm

Purification efficiency: 100(U.S. Federal Standard 209E)

Weight: 85kg

The average number of colonies: ?0.5/dish, when

Noise: ?62db(A)

Half peak vibration: ?3?m(x*y*z)

Power: 220V50HZ

Power consumption: 0.3kw

The average wind speed: 0.3-0.6 m/s(adjustable)

Fitter size: 695*455*50*1(L*W*H) mm

Lighting/germicidal lamp: 15w*1*15w*1

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,