

Product Name :
Urine Analyzer System Test Machine

Product Code :
CLNCL-CDD0005



Description :

Urine Analyzer System Test Machine

Technical Specification :

Test parameters: URS-10T, URS-11
wavelength: 558nm, 635nm, 720nm
Test speed: Up to 120 tests/hour at Continuous
Memory Size: 1000 test results
optional : 10 or 11 parameters
Print: Thermal printer inside, with outer stylus printer interface
Host-communication: PC connected with RS232
Environment Condition: 0-40°C, RH Power Supply: AC 220V(±15%), 50-60Hz
Power: ?60W