

Product Name :
Laboratory Oven (Hot Air Oven)**Product Code :**
NMCBS11040**Description :**

Laboratory Oven (Hot Air Oven) We offer laboratory ovens that are used in a variety of thermal processing applications including general lab work, component and stability testing, core hardening, drying glassware, and sterilizing. These are known for their robust construction, low maintenance and high efficiency. Some of the features are as follows: •Double Walled Construction •Exterior & Interior made of stainless steel 304/316 (GMP model) or Exterior made of Mild Steel with powder coating finish and interior made of stainless steel 304 (STD Model) •80/20 Ni Heating Elements placed at bottom •Temp. Range : 5°C above ambient to 300°C 2°C •230V AC 50HZ Mains Supply. •Digital Temperature Controller. •Fan arrangement for forced air circulation. Models Chamber H x W x D (cms) Watts Shelves AEI-20a 30 x 30 x 30 750 1 AEI-20b 35 x 35 x 35 1000 2 AEI-20c 45 x 45 x 45 1500 2 AEI-20d 45 x 45 x 60 1750 2 AEI-20e 60 x 60 x 60 2000 3 AEI-20f 45 x 60 x 90 2200 3 AEI-20g 60 x 60 x 90 2500 3 Laboratory Oven (Hot Air Oven)