

Product Name :
Pathology Lab Semi Auto Analyzer**Product Code :**
MED-PTH0016**Description :**

Pathology Lab Semi Auto Analyzer

Technical Specification :

The Semi-Auto Analyzer helps in measuring chemical indexes by examining blood and other body fluid. Further, it is combined with other clinical information, to assist in evaluating organs function, diagnosing diseases, recognizing disease gene and determining the norm for future therapy. The in-built thermal printer of the analyzer prints the results and makes, diagnosing a disease simpler process.

Product Details :

Application: Laboratory Use
Display Type: TFT touch-screen
Display Size: 7.0 inch
Number of Test Modes: 2
Test Mode Type: Flowcell or cuvette
Packaging Type: Boxes

SALIENT FEATURES

7.0" TFT touch-screen & pop-up keypad
Two test mode: flowcell or cuvette
Supports bi-chromatic tests for end point, fixed-time, kinetics methods
External keyboard or mouse support via USB
User-friendly interface and multi-language operation software
Built-in thermal printer

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE,](#)